

# PaxeraVet

Veterinary image acquisition and viewing station

PaxeraVet is a full-featured browser-based PACS solution to acquire, view, process, distribute and manage veterinary DICOM data.

The thin-client multi-modality medical viewer features comprehensive diagnostics and reporting tools that provides powerful imaging functions necessary to today's veterinary practices. The solution is very flexible and easy to manage, its intuitive interface helps users to seamlessly use the application with less training and technical knowhow.



## Key features include:-

- Supports all DICOM modalities,
- DICOM import/export function of any media
- Work list function to register and retrieve animal name and owner demographics information.
- Work list with veterinary specific information including breed and studbook names
- User-friendly interface designed to provide easy and streamlined access to images and reports
- Image processing tools and filters including invert, rotate, flip, resize, emboss and median
- Comprehensive measurement and annotation tools
- Region of interest (ROI) processing and marking
- Easy to use filming tool, print selected images on DICOM or Windows printers
- Customizable toolbar for the most frequent used functions
- Archive studies to CD/DVD, with an embedded viewer for viewing studies and reports anywhere
- Powerful image display with advanced hanging protocols and individual preferred viewing interface
- Extensive security features including user access control and automatic browser logout
- Zero footprint viewer enabling veterinarians to view studies from any browser or mobile device
- Concurrently open different studies for comparison
- Reporting and transcription tool with image insertion and auto header
- Optional cloud hosting and backup for a nominal fee-per-study
- iPad/ iPhone access with report dictation via iPaxera App
- Interfaces with Veterinary Information systems via PaxeraBroker®
- Flexible Licensing options for 2, 5, 10 and unlimited modalities